

Human CITED2 Recombinant Protein

Catalog Number: PROTQ99967

Overview

Product Name	Human CITED2 Recombinant Protein
Description	Human CITED2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-270aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 50% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q99967, 1-270aa; MADHMMAMNH GRFPDGTNGL HHHPAHRMGM GQFSPHHHQ QQPQHAFNA LMGEHIHYGA GNMNATSGIR
HAMGPGTVNG GHPPSALAPA ARFNNSQFMG PPVASQGGS L PASMQLQKLN NQYFNHHPYP HNHYMPDLHP AAGHQMNGTN QHFRDCNPKH
SGGSSTPGGS GGSSTPGGS SSSGGGAGSS NSGGGSGSN MPASVAHVPA AMLPPNVIDT DFIDEEV LMS LVIEMGLDRI KELPELWLGQ
NEFDFTDFV CKQQPSRVSC LEHHHHHH [?/description_after_attributes?](#)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CITED2 Recombinant Protein