

Human CITED2 Recombinant Protein

Catalog Number: PROTQ99967

Overview

| | |
|--------------|--|
| Product Name | Human CITED2 Recombinant Protein |
| Description | Human CITED2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-270aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size | Starting from 20ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 50% glycerol, 1mM DTT |

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q99967, 1-270aa; MADHMMAMNH GRFPDGTNGL HHHPAHRMGM GQFSPHHHQ QQPQHAFNA LMGEHIHYGA GNMNATSGIR
HAMGPGTVNG GHPPSALAPA ARFNNSQFMG PPVASQGGSL PASMQLQKLN NQYFNHHYPY HNHYPDLHP AAGHQMNGTN QHFRDCNPKH
SGGSSTPGGS GGSSTPGGSG SSSGGGAGSS NSGGGSGSGN MPASVAHVPA AMLPPNVIDT DFIDEEVLMS LVIEMGLDRI KELPELWLGQ
NEFDFTDFV CKQQPSRVSC LEHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CITED2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.