

## Mouse PBEF/Visfatin Recombinant Protein

Catalog Number: PROTQ99KQ4

### Overview

Product Name	Mouse PBEF/Visfatin Recombinant Protein
Description	Mouse PBEF/Visfatin Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-491aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4)

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q99KQ4, 1-491aa; MNAAAEAEFN ILLATDSYKV THYKQYPPNT SKVYSYFECR EKTENSKVR KVKYEETVFY GLQYILNKYL KGKVVTKEDI  
QEAKVYREH FQDDVFNERG WNYILEKYDG HLPYEVKAVP EGSVIPRGNV LFTVENTDPE CYWLTNWIET ILVQSWYPIT VATNSREQKK  
ILAKYLLETS GNLDGLEYKL HDFGYRGVSS QETAGIGASA HLVNFKGTDV VAGIALIKKY YGTDKDPVPGY SVPAAEHSTI TAWGKDHEKD  
AFEHIVTQFS SVPVSVSDS YDIYNACEKI WGEDLRHLIV SRSTEAPLII RPDGSGNPLDT VLKVLIDILGK KFPVTENSKG YKLLPPYLRV IQGDGVDINT  
LQEIVEGMKQ KKWSIENVSF GSGGALLQKL TRDLLNCSFK CSYVVTNGLG VNVFKDPVAD PNKRSKKGRLL SLHRTPAGNF VTLEEGKDDL  
EEYGHLLHT VFKNGKVTKS YSFDEVRKNA QLNIEQDVAP H

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse PBEF/Visfatin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.