

Human 17 beta -HSD14/HSD17B14 Recombinant Protein

Catalog Number: PROTQ9BPX1

Overview

Product Name	Human 17 beta -HSD14/HSD17B14 Recombinant Protein
Description	Human 17 beta -HSD14/HSD17B14 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-270aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q9BPX1, 1-270aa; MATG TRYAGKVVVV TGGGRGIGAG IVRAFNNSGA RVVICDKDES GGRALEQELP GAVFILCDVT QEDDVKTLVS ETIRRFGRDL CVVNNAGHHP PPQRPEETSA QGFRQLELN LLGTYTLTKL ALPYLRKSQG NVINISLVG AIGQAQAVPY VATKGAVTAM TKALALDESP YGVRVNCISP GNIWTPWEE LAALMPDRA TIREGMLAQP LGRMGQPAEV GAAAVFLASE ANFCTGIELL VTGGAELGYG CKASRSTPVD APDIPS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human 17 beta -HSD14/HSD17B14 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.