

Human DPH2 Recombinant Protein

Catalog Number: PROTQ9BQC3

Overview

Product Name	Human DPH2 Recombinant Protein
Description	Human DPH2 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 1-489aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 1mM DTT, 20% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q9BQC3, 1-489aa; MESMFSSPAE AALQRETGVP GLLTLPDLD GVEYELERVAG FVRDLGCERV ALQFPDQLLG DAVAVAARLE
ETTGSKMFI L GDTAYGSCCV DVLGAEQAGA QALIHFGPAC LSPPARPLPV AFVLRQRSVA LELCVKAFEQA QNPDPKAPVV LLSEPACAH
LEALATLLRP RYLDLLVSSP AFPQPVGSL S PEPMLERFG RRFPLAPGRR LEEYGAFYVG GSKASDPDL DPDLRLLLG WAPGQPFSSC
CPDTGKTQDE GARAGRLRAR RRYLVERARD ARVVGLLAGT LGVAQHREAL AHLRNLTQAA GKRSYVLALG RPTPAKLANF PEVDVFLLA
CPLGALAPQL SGFFQPILA PCELEAACNP AWPPPGLAPH LTHYADLLPG SPFHVALPPP ESELWETPDV SLITGDLRPP PAWKSSNDHG
SLALTPRPQL ELAESSPAAS FLSSRSWQGL EPRLGQTPVT EAVSGRRGIA IAYEDEGSG

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human DPH2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.