

## Human TIPIN Recombinant Protein

Catalog Number: PROTQ9BVW5

### Overview

Product Name	Human TIPIN Recombinant Protein
Description	Human TIPIN Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-301aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q9BVW5, 1-301aa; MLEPQEN GVIDLPDYEH VEDEFPPFP PPASPERQDG EGTEPDEESG NGAPVPVPPK RTVKRNIPKL DAQLISERG LPALRHVFDK AKFKGKGHEA EDLKMLIRHM EHWAHRLF PK LQFEDFIDRV EYLGSKKEVQ TCLKRIRLDL PILHEDFVSN NDEVAENNEH DVTSTELDPF LTNLSESEMF ASELRSRLTE EQQQRIERNK QLALERRQAK LLSNSQTLGN DMLMNTPRAH TVEEVNTDED QKEESNGLNE DILDNPCNDA IANTLNEEET LLDQSFKNVQ QQLDATSRNI TEAR

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human TIPIN Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.