

Human recombinant CTRP5/C1qTNF5 protein

Catalog Number: PROTQ9BXJ0-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant CTRP5/C1qTNF5 protein
Description	Human recombinant CTRP5/C1qTNF5 protein expressed in E.coli with His-Tag. Sequence domain: 16-243 aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E.coli
Formulation	20mM Tris-HCl buffer (pH8.5) containing 30% glycerol, 0.2M NaCl

Concentration

0.25mg/ml (determined by BCA assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q9BXJ0, 16-243 aa; MGSSHHHHHH SGLVPRGSH MGSMSPPLD DNKIPSLCPG HPGLPGTPGH HGSQGLPGRD GRDGRDGAPG
APGEKGEGGR PGLPGPRGDP GPRGEAGPAG PTGPAGECSV PPRSFAFAKR SESRVPPPSD APLPFDRVLV NEQGHYDAVT GKFTCQVPGV
YYFAVHATVY RASLQFDLVK NGESIASFFQ FFGGWPKPAS LSGGAMVRLE PEDQVWVQVG VGDIYIGIYAS IKTDSTFSGF LVYSDWHSSP VFA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant CTRP5/C1qTNF5 protein](#)

For Research Use Only. Not for use in diagnostic procedures.