

## Human PECR Recombinant Protein

Catalog Number: PROTQ9BY49

### Overview

|              |   |
|--------------|---|
| Product Name | Human PECR Recombinant Protein  |
| Description  | Human PECR Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-303aa. Application(s): SDS-PAGE. |
| Size         | Starting from 20ug  |
| Tag          | His-Tag   |
| Form         | Liquid  |
| Source       | E. coli   |
| Formulation  | 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2mM DTT, 0.15M NaCl.   |

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, Q9BY49, 1-303aa; MGSSHHHHHH SSSLVPRGSH MGSMMASWAK GRSYLAPGLL QGQVAIVTGG ATGIGKAIVK ELLELGSNVV IASRKLERLK SAADELQANL PPTKQARVIP IQCNIRNEEE VNNLVKSTLD TFGKINFLVN NGGGQFLSPA EHISSKGWHA VLETNLTGTF YMCKAVYSSW MKEHGGSIVN IIVPTKAGFP LAVHSGAARA GYYNLTKSLA LEWACSGIRI NCVAPGVIYS QTAVENYGSW GQSFEGSFQ KIPAKRIGVP EEVSSVVCFL LSPAASFITG QSVDDVDGGRS LYTHSYEVPD HDNWPKGAGD LSVVKKMKET FKEKAKL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PECR Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.