

Human IGSF4A/SynCAM1 Recombinant Protein

Catalog Number: PROTQ9BY67

Overview

Product Name	Human IGSF4A/SynCAM1 Recombinant Protein
Description	Human IGSF4A/SynCAM1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 45-374aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q9BY67, 45-374aa; MGSSHHHHHH SSSLVPRGSH MGSQNLFTKD VTVIEGEVAT ISCQVNKSDDD SVIQLLNPNR QTIIYFRDFRP
LKDSRFQLLN FSSSELKVSL TNVSISDEGR YFCQLYDPP QESYTTITVL VPPRNLMIDI QKDTAVEGEE IEVNCTAMAS KPATTIRWFK GNTLKGKSE
VEEWSDMYTV TSQMLMLKVHK EDDGVPVICQ VEHPAVTGNL QTQRYLEVQY KPQVHIQMTY PLQGLTREGD ALELTCEAIG KPQPVMVTWV
RVDDEMPQHA VLSGPNLFIN NLNKTDNGTY RCEASNIVGK AHSDYMLVYV DPPTTIPPT TTTTTTTTTT TTILTIITDS RAGEEGSIRA VDH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human IGSF4A/SynCAM1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.