

Human recombinant RAB34 protein

Catalog Number: PROTQ9BZG1-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant RAB34 protein
Description	Human recombinant RAB34 protein expressed in E.coli with His-Tag. Sequence domain: 1-259aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E.coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 40% glycerol, 0.1M NaCl, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9BZG1, 1-259aa; MGSSHHHHHH SGLVPRGSH MNILAPVRRD RVLAE LPQCL RKEAALHG HK DFHPRVTCAC QE HRTGTVGF
KISKVIVVGD LSVGKTCLIN RFCKDTFDKN YKATIGVDFE MERFEVLGIP FSLQLWDTAG QERFKCIAST YYRGAQAI III VFNLNDVASL
EHTKQWLADA LKENDPSSVL LFLVGSKKDL STPAQYALME KDALQVAQEM KAEYWAVSSL TGENVREFF F RVAALTFEAN VLAELEKSGA
RRIGDVVRIN SDDSNLYLTA SKKKPTCCP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant RAB34 protein](#)

For Research Use Only. Not for use in diagnostic procedures.