

Human uCK2 Recombinant Protein

Catalog Number: PROTQ9BZX2

Overview

| | |
|--------------|---|
| Product Name | Human uCK2 Recombinant Protein |
| Description | Human uCK2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-261aa. Application(s): SDS-PAGE. |
| Size | Starting from 20ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 2mM DTT, 100mM NaCl, 0.1mM PMSF, 1mM EDTA |

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9BZX2, 1-261aa; MAGDSEQLQ NHQQPNGGEP FLIGVSGGTA SGKSSVCAKI VQLLGQNEVD YRQKQVVILS QDSFYRVLTS
EQKAKALKGQ FNFDPDAFD NELILKTLKE ITEGKTQVIP VYDFVSHSRK EETVTVPAD VVLFEGILAF YSQEVRLDFQ MKLFDVTDAD
TRLRRVLRD ISERGRDLEQ ILSQYITFVK PAFEEFCLPT KKYADVIIPR GADNLVAINL IVQHIQDILN GGPSKRQTNG CLNGYTPSRK RQASESSRP
H

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human uCK2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.