

TNKS1BP1 (NM_033396) Human Recombinant Protein

Catalog Number: PROTQ9C0C2

Overview

Product Name	TNKS1BP1 (NM_033396) Human Recombinant Protein
Description	Recombinant protein of human tankyrase 1 binding protein 1, 182kDa (TNKS1BP1)
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at -80°C.

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



TNKS1BP1 (NM_033396) Human Recombinant Protein