

## Mouse Cyclophilin-40/Ppid Recombinant Protein

Catalog Number: PROTQ9CR16

### Overview

Product Name	Mouse Cyclophilin-40/Ppid Recombinant Protein
Description	Mouse Cyclophilin-40/Ppid Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-370aa. Application(s): SDS-PAGE, Enzyme Activity.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT

### Concentration

1mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q9CR16, 1-370aa; MSHAS PAAKPSNSKN PRVFFDVDIG GERVGRIVLE LFADIVPKTA ENFRALCTGE KGTGSTTGKP LHFKGCPFHR IIKKFMIQGG DFSNQNGTGG ESIYGEKFED ENFHVKHRE GLLSMANAGP NTNGSQFFIT TVPTPHLDGK HVVFGQVIKG LGVARTLENV EVNGEKPAKL CVIAECGELK EGDDWGIFPK DSGSDSHPDF PEDADIDLKD VDKILLISED LKNIGNTFFK SQNWEMAIAKK YAKVLRVDS SKAVIEKADR SRLQPIALSC VLNIGACKLK MSNWQGAIDS CLEALEMDPS NTKALYRKAQ GWQGLKEYDQ ALADLKKAAQE IAPGDKAIQA ELLKVKQMIK AQKDKKAVY AKMFA

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Cyclophilin-40/Ppid Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.