

Mouse PGAM1 Recombinant Protein

Catalog Number: PROTQ9DBJ1

Overview

Product Name	Mouse PGAM1 Recombinant Protein
Description	Mouse PGAM1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-254aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, Q9DBJ1, 1-254aa; MGSSHHHHHH SGLVPRGSH MGSMAAYKL VLIRHGSAW NLENRFSGWY DADLSPAGHE EAKRGGQALR
DAGYEFDICF TSVQKRAIRT LWTVLDAIDQ MWLPVVRTWR LNERHYGGLT GLNKAETA AK HGEA QVKIWR RSYDVPPPPM EPDHPFYSNI
SKDRRYADLT EDQLPSCESL KDTIARALPF WNEEIVPQIK EGKRVLIAAH GNSLRGIVKH LEGLSEEAIM ELNLPTGIPI VYELDKNLKP IKPMQFLGDE
ETVRKAMEAV AAQGKVKK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse PGAM1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.