

## Mouse MESDC2 Recombinant Protein

Catalog Number: PROTQ9ERE7

### Overview

Product Name	Mouse MESDC2 Recombinant Protein
Description	Mouse MESDC2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 30-224aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

1mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q9ERE7, 30-224aa; MGSSHHHHHH SGLVPRGSH MGSADTPGEA TPPPRKKKDI RDYNDADMAR LLEQWEKDDD IEEGDLPEHK  
RPSAPIDFSK LDPGKPESIL KMTKKGKTL M FVTVSGNPT EKETEEITSL WQGS LFNANY DVQRFIVGSD RAIFMLRDGS YAW EIKDFLV  
SQDRCAEVT L EGQMYPGKGG GSKEKNKTKP EKAKKKEGDP KPRASKEDNR AGSRREDL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse MESDC2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.