

Human SRR Recombinant Protein

Catalog Number: PROTQ9GZT4

Overview

Product Name	Human SRR Recombinant Protein
Description	Human SRR Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-340aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q9GZT4, 1-340aa; MGSSHHHHHH SGLVPRGSH MGS HMCAQYC ISFADVEKAH INIRDSIHLT PVL TSSILNQ LTGRNLFFKC
ELFQKTGSFK IRGALNAVRS LVPDALERKP KAVVTHSSGN HGQALTYAAK LEGIPAYIVV PQTAPDCKKL AIQAYGASIV YCEPSDESRE
NVAKRVT EET EGIMVHPNQE PAVIAGQGTI ALEVLNQVPL VDALVVPVGG GGMLAGIAIT VKALKPSVKV YAAEPSNADD CYQSKLKGKL
MPNLYPPE TI ADGVKSSIGL NTWPIIRDLV DDIFTVTEDE IKCATQLVWE RMKLLIEPTA GVGVAAVLSQ HFQTVSPEVK NICIVLSGGN
VDLTSSITWV KQAERPASYQ SVSV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SRR Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.