

## Human TAZ/WWTR1 Recombinant Protein

Catalog Number: PROTQ9GZV5

### Overview

Product Name	Human TAZ/WWTR1 Recombinant Protein
Description	Human TAZ/WWTR1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-400aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.4M uREA, 10% glycerol

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q9GZV5, 1-400aa; MNPASAP PPLPPPGQQV IHVTQDLDTD LEALFNSVMN PKPSSWRKKI LPESFFKEPD SGSHSRQSST DSSGGHPGPR LAGGAQHVRS HSSPASLQLG TGAGAAGSPA QQHAHLRQQS YDVTDELPLP PGWEMTFTAT GQRYFLNHIE KITTWQDPRK AMNQPLNHMN LHPAVSSTPV PQRSMASVQP NLVMNHQHQQ QMAPSTLSQQ NHPTQNPPAG LMSMPNALT TQQQQQKLRL QRIQMERERI RMRQEELMRQ EAALCRQLPM EAETLAPVQA AVNPPTMPD MRSITNNSD PFLNGGPYHS REQSTD SGLG LGCYSVPTTP EDFLSNVDEM DTGENAGQTP MNINPQQTRF PDFLDCLPGT NVDLGTLESE DLIPLFNDVE SALNKSEPFLL TWL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human TAZ/WWTR1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.