

Human CDCP1 Recombinant Protein

Catalog Number: PROTQ9H5V8

Overview

Product Name	Human CDCP1 Recombinant Protein
Description	Human CDCP1 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 30-667aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9H5V8, 30-667aa; FEIALPRESN ITVLIKLGTP TLLAKPCYIV ISKRHITMLS IKSGERIVFT FSCQSPENHF VIEIQKNIDC MSGPCPFGEV
QLQPSTSLLP TLNRTFIWDV KAHKSIGLEL QFSIPRLRQI GPGESC PDGV THSISGRIDA TVVRIGTFCS NGTVSRIKMQ EGVKMALHLP
WFHPRNVSGF SIANRSSIKR LCIIESVFEG EGSATLMSAN YPEGFPEDEL MTWQFVVP AH LRASVSFLNF NLSNCERKEE RVEYYIPGST
TNPEVFKLED KQPGNMAGNF NLSLQGCDDQ AQSPGILRLQ FQVLVQHPQN ESNKIYVVDL SNERAMSLTI EPRPVKQSRK FVPGCFVCLC
SRTCSSNLTL TSGSKHKISF LCDDLTRLWM NVEKTISCTD HRYCQRKSYS LQVPSDILHL PVELHDFSWK LLVPKDRLSL VLVPAQKLQ
HTHEKPCNTS FSYLVASAIP SQDLYFGSFC PGGSIKIQV QKNISVTLRT FAPSFQQEAS RQGLTVSFIP YFKEEGVFTV TPDTKSKVYL
RTPNWDRGLP SLTSVSWNIS VPRDQVACLT FFKERSGVVC QTGRAFMIIQ EQRTRAEEIF SLDEDVLPKP SFHHHSFWVN ISNCSPTSGK
QLDLLFSVTL TPRTVDLT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CDCP1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.