

C19orf44 (NM_032207) Human Recombinant Protein

Catalog Number: PROTQ9H6X5

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	C19orf44 (NM_032207) Human Recombinant Protein
Description	Recombinant protein of human chromosome 19 open reading frame 44 (C19orf44)
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 µg/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MASARKASRPMRDVFGDFSDVSLEDSTMEEIRNFQISRNLTKIAPGHSRFLKRNQTLDEKHLKKENPVLGSGPRLASCRPPTTASRIRANAALMKLAQLE
TRIMNRKLQRNLSDTESDSMTADAGLPKRADRLSGGALELASQNTDKTSQNTQARELPVTENNAQNAKVSRLFKKKQAPVENISPEAPAGKERTLQTPK
QKEPARTFDSPDSDEEEMKVLLGSLMDSSREKNTNQGFSSANVSEEEERKLFSVPSQLRAFTVPVELSSAKPSQTSHLPTSLAADRTLHSTRSRADYPQS
HVSSDTASHTPSVSITGAFSNSVSLKMGHVKLVSPPGRSEAEVDEPVSEGADDSLDEFNRINILSLDGLAPAVSENSDLEQEEESAQRQKTAGKIFRAEAST
GQDAPRQAQARSWASQGAASAEGDESEVSEHLSASSASAIQQDSTSSMQPPSEAPMVNTVSSAYSEDFENSPSLTASEPTAHSKESLDRTLDALSESSSS
VKTDLPQTAESRKKSGRHVTRVLVKDTAVQTPDPAFTYEWTKVASMAAMGPALGGAYVDPTPIANHVISADAIEALTAYSPAVLALHDVLKQQLSLTQ

QFIQASRHLHASLLRSLDADSFHYHTLEEAKEYIRCHRPAPLTMEDALEEVNKL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

C19orf44 (NM_032207) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.