

## SYNPO2L (NM\_001114133) Human Recombinant Protein

Catalog Number: PROTQ9H987

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|              |                                                                                                  |
|--------------|--------------------------------------------------------------------------------------------------|
| Product Name | SYNPO2L (NM_001114133) Human Recombinant Protein                                                 |
| Description  | Purified recombinant protein of Homo sapiens synaptopodin 2-like (SYNPO2L), transcript variant 1 |
| Size         | 20 µg                                                                                            |
| Tag          | C-Myc/DDK                                                                                        |
| Form         | Frozen Solution in PBS Buffer                                                                    |
| Source       | HEK293T                                                                                          |
| Formulation  | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol                                             |

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MGAEEEVLVTLSGGAPWGFRLHGGAEQRKPLQVSKIRRRSQAGRAGLRERDQLLAINGVSCTNLSHASAMSLIDASGNQLVLTQRLADEGPVQSPSPH  
ELQVLSPLSPLSPEPPGAPVPQPLQPGSLRSPDSEAYYGETDSDADGPATQEKPRRPRRRGPTRPTPPGAPPDEVYLSDDSPAEPAPTIPGPPSQGDSRV  
SPSWEDGAALQPPPAEALLPHGPLRPGPHLIPMVGVPVPHVAEDLTITYTQKAKQAKLQRAESLQEKSIAEAKTKCRTIASLLTAAPNPHSKGVL MFKK  
RRQRAKKYTLVSFGAAAGTGAEEDGVPPTSESELDEEAFSDARSLTNQSDWDSPLYDMELARAGSRASEGQGSGLGGQLSEVSGRGVQLFEQQRQRA  
DSSTQELARVEPAAMLNGEGLQSPPRASAPPEAAVLPPSPLPAPVASPRPFQPGGGAPTAPSFNRSARPTPLGQGRPTTTSVIFRPLAPKRANDSL  
GGLSPAPPPFLSSQGPTPLSFTSGVPSHAPVSGSPSTPRSSGPVTATSSLYIPASRPVTPGGAPEPPAPPSAAAMTSTASIFLSAPLRPSARPEAPAPGPGA

PEPPSAREQRISVPAARTGILQEARRRGTRKQMFRPGKEETKNSPNELLSLVQNLDEKPRAGGAESGPEEDALSLGAEACNFMQPVGARSYKTLPHVTP  
KTPPPMAPKTPPPMTPKTPPPVAPKPPSRGLLDGLVNGAASSAGIPEPPRLQGRGGELFAKRQSRADRYVVEGTPGPGLGPRPRSPSPTPSLPPSWKYSP  
NIRAPPPAIYNPLLSPPFPQAARTLPKAQSQGPRAATPKQGKALDFMRHQPYQLKTAMFCFDEVPTPGPIASGSPKTARVQEIRRFSTPAPQPTAEPLAP  
TVLAPRAATTLDDEPIWRTELASAPVPSAPPPEAPRGLGASPSGCFQVARPRFSATRTGLQAHVWRPGAGHQ

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[SYNPO2L \(NM\\_001114133\) Human Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.