

Human PPCS Recombinant Protein

Catalog Number: PROTQ9HAB8

Overview

Product Name	Human PPCS Recombinant Protein
Description	Human PPCS Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-311aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q9HAB8, 1-311aa; MGSSHHHHHH SGLVPRGSH MAEMDPVAEF PQPPGAARWA EVMARFAARL GAQGRRVVLV TSGGTKVPLE
ARPVRFIDNF SSGRRGATSA EAFLAAGYGV LFLYRARSF PYAHRFPPT WLSALRPSGP ALSGLLSLEA EENALPGFAE ALRSYQEA
AGTFLAVEFT TLADYLHLLQ AAAQALNPLG PSAMFYLA VSDFYVPVSE MPEHKIQSSG GPLQITMKMV PKLLSPLVKD WAPKAFIISF
KLETDP AIVI NRARKALEIY QHQVVANIL ESRQSFVIV TKDSETKLLL SEEEIEKGV EIEKIVDNLQ SRHTAFIGDR N

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PPCS Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.