

EML4 (NM_019063) Human Recombinant Protein

Catalog Number: PROTQ9HC35

Overview

Product Name	EML4 (NM_019063) Human Recombinant Protein
Description	Purified recombinant protein of Homo sapiens echinoderm microtubule associated protein like 4 (EML4), transcript variant 1
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 µg/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MDGFAGSLDDSSISAASSTSDVQDRLSALESRVQQQEDEITVLKAAALADVLRRLAISEDHVASVKKSVSSKGQPSRAVIPMSCITNGSGANRKPSTSAVSI
AGKETLSSAAKSGTEKKKEKPQGQREKKEESHNSDQSPQIRASPSQPSSQPLQIHRQTPESKNATPTKSIKRPSPAEKSHNSWENSDDSRNKLKIPSTP
KLIPKVTKTADKHKDVIINQEGEYIKMFMGRGPITMFI PSDVDNYDDIRTEL PPEKLEWAYGYRGKDCRANVYLLPTGEIVYFIASVVVLFNYEERTQR
HYLGHTDCVKCLAIHPDKIRIATGQIAGVDKDGRLQPHVRVWDSVTLSTLQIIGLGTFERGVGCLDFSKADSGVHLCVIDDSNEHMLTVWDWQKKA
GAEIKTTNEVVLAVEFHPTDANTIITCGKSHIFFWTWSGNSLTRKQGIFGKYEKPKFVQCFLAFLNGDVLTDGSDGGVMLIWSKTTVEPTPGKGPVYQI
SKQIKAHDGVSFTLCQMRNGMLLTGGGKDRKIILWDHDLNPEREIEVPDQYGTIRAVAEGKADQFLVGTSRNFILRGTFNDGFQIEVQGHTDELWGLA
THPFKDLLLTCADRQVCLWNSMEHRLEWTRLVDEPGHCADFHPSGTVAIGTHSGRWFVLD AETRDLVSIHTDGNEQLSVMRYSIDGTFLAVGSHD
NFIYLYVVENGRKYSRYGRCTGHSSYITHLDWSPDNKYIMNSGDYEILYWDIPNGCKLIRNRSDCKDIDWTTYTCVLGFQVFGVWPEGSDGTDINAL
VRSHNRKVIAVADDFCKVHLFQYPCSKAKAPSHKYSAHSSHVTVNSFTHNSHLISTGGKDMSIQWKLVEKLSLPQNETVADTTLT KAPVSSTESVIQS
NTPTPPSQPLNETAEESRISSTLLENSLEQTVEPSEDHSEEESEEGSGDLGEPLEYEPCNEISKEQAKATLLEDQQDPSPSS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



EML4 (NM_019063) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.