

Mouse IL-17 RB Recombinant Protein

Catalog Number: PROTQ9JIP3

Overview

Product Name	Mouse IL-17 RB Recombinant Protein
Description	Mouse IL-17 RB Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 18-286aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, Q9JIP3, 18-286aa; ADPREPTIQC GSETGPSPEW MVQHTLTPGD LRDLQVELVK TSVAEEFSI LMNISWILRA DASIRLLKAT KICVSGKNNM NSYSCVRCNY TEAFQSQTRP SGGKWTFSYV GFPVELSTLY LISAHNIPNA NMNEDSPSLS VNFTSPGCLN HVMKYKKQCT EAGSLWDPDI TACKKNEKMOV EVNFTTNPLG NRYTILIQRD TTLGFSRVLE NKLMRTSVAI PVTEESEGAV VQLTPYLHTC GNDCIRREGT VVLCSETSAP IPPDDNRRML GGVEPKSCDK THTCPPCPAP ELLGGPSVFL FPPKPKDTLM ISRTPEVTCV VVDVSHEDPE VKFNWYVDGV EVHNAKTKPR EEQYNSTYRV VSVLTVLHQD WLNGKEYKCK VSNKALPAPI EKTISKAKGQ PREPQVYTLPSRDELTKNQ VSLTCLVKGF YPSDIAVEWE SNGQPENNYK TTPPVLDSDG SFFLYSKLTV DKSRWQQGNV FSCVMHEAL HNHYTQKSL SLPQKHHHHH H

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-17 RB Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.