

Mouse recombinant BACE-2 protein

Catalog Number: PROTQ9JL18

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse recombinant BACE-2 protein
Description	Mouse recombinant BACE-2 protein expressed in HEK293 with His-Tag. Sequence domain: 20-462aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, Q9JL18, 20-462aa; AVPALAPAPF TLPLQVARAT NHRASAVPGL GTPELPRADG LALALEPVRA TANFLAMVDN LQGDSGRGYY
LEMLIGTPPQ KVQILVDTGS SNFAVAGAPH SYIDTYFDSE SSSTYHSKGF DVTVKYTQGS WTGFVGEDLV TIPKGFNSSF LVNIATIFES
ENFFLPGIKW NGILGLAYAA LAKPSSSLET FFDLSVAQAK IPDIFSMQMC GAGLPVAGSG TNGGSLVLGG IEPsLYKGDI WYTPIKEEWY YQIEILKLEI
GGQNLNLCR EYNADKAIVD SGTTLRLRPQ KVFDAAVVEAV ARTSLIPEFS DGFWTGAQLA CWTNSETPWA YFPKISIIYLR DENASRSFRI

TILPQLYIQP MMGAGFNyec YRFGISSSTN ALVIGATVME GFYVVFDRaq RRVGFAVSPC AEIEGTTVSE ISGPFSTEDI ASNCVPAQAL NEP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Mouse recombinant BACE-2 protein

For Research Use Only. Not for use in diagnostic procedures.