

## Mouse IL-36 alpha/IL-1F6 Recombinant Protein

Catalog Number: PROTQ9JLA2

### Overview

Product Name	Mouse IL-36 alpha/IL-1F6 Recombinant Protein
Description	Mouse IL-36 alpha/IL-1F6 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-160aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

1mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q9JLA2, 1-160aa; MNKEKEL RAASPSLRHV QDLSSRVWIL QNNILTAVPR KEQTVPTIT LLPCQYLDLTL ETNRGDPTYM GVQRPMSCLF CTKDGEQPV LQLGEGNIMEM YNKKEPVKAS LFYHKKSGTT STFESAAPFG WFIIVCSKGS CPLILTQELG EIFITDFEMI VVH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse IL-36 alpha/IL-1F6 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.