

Human recombinant DC-SIGN/CD209 protein

Catalog Number: PROTQ9NNX6

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant DC-SIGN/CD209 protein
Description	Human recombinant DC-SIGN/CD209 protein expressed in Baculovirus with hIgG-His-Tag. Sequence domain: 60-404aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9NNX6, 60-404aa; VSKVPSS ISQEQSRQDA IYQNLTLQKA AVGELSEKSK LQEIQELTQ LKAAVGELPE KSKLQEIQE LTRLKAAVGE LPEKSKLQEI YQELTWLCAA VGELPEKSKM QEIQELTRL KAAVGELPEK SKQQEIQEL TRLKAAVGEL PEKSKQQEIQ QELTRLKAAV GELPEKSKQQ EIQELTQLK AAVERLCHPC PWEWTFQGN CYFMSNSQRN WHDSITACKE VGAQLVVIKS AEEQNFLQLQ SSRSNRFTWM GLSDLNQEGT WQWVDGSPLL PSFKQYWNRG EPNNVGEEDC AEFSGNGWND DKCNLAKFWI CKKSAASCSR DEEQFLSPAP ATPNPPPA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human recombinant DC-SIGN/CD209 protein

For Research Use Only. Not for use in diagnostic procedures.