

## Human FBXO6 Recombinant Protein

Catalog Number: PROTQ9NRD1

### Overview

Product Name	Human FBXO6 Recombinant Protein
Description	Human FBXO6 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-293aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q9NRD1, 1-293aa; MGSSHHHHHH SGLVPRGSH MGSMDAPHSK AALDSINELP ENILLELFTH VPARQLLLNC RLVCSLWRDL  
IDLMTLWKRK CLREGFITKD WQPVADWKI FYFLRSLHRN LLRNPCAED MFAWQIDFNG GDRWKVESLP GAHGTFDPDP KVKKYFVTSY  
EMCLKSQLVD LVAEGYWEEL LDTRPDIVV KDWFAARADC GCTYQLKVQL ASADYFVLAS FEPPPVTIQQ WNNATWTEVS YTFSDYPRGV  
RYILFQHGGGR DTQYWAGWYG PRVTNSSIVV SPKMTRNQAS SEAQPGQKHG QEAAQSPYR AVVQIF

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human FBXO6 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.