

Human CCL28 Recombinant Protein

Catalog Number: PROTQ9NRJ3

Overview

Product Name	Human CCL28 Recombinant Protein
Description	Human CCL28 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 23-127aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	10mM Sodium Citrate buffer (pH 3.5) containing 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9NRJ3, 23-127aa; MGSSHHHHHH SGLVPRGSH MILPIASSCC TEVSHHISRR LLERVNMCRI QRADGDCDLA AVILHVKRRR
ICVSPHNHTV KQWMKVQAAK KNGKGNVCHR KKHGKRNSN RAHQGKHETY GHKTPY [?/description_after_attributes?](#)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CCL28 Recombinant Protein