

## Human RASSF1A Recombinant Protein

Catalog Number: PROTQ9NS23

### Overview

Product Name	Human RASSF1A Recombinant Protein
Description	Human RASSF1A Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-340aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q9NS23, 1-340aa; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMSE PELIELRELA PAGRAGKGRT RLERANALRI  
ARGTACNPTR QLVPGRGHRF QPAGPATHTW CDLCGDFIWG VVRKGLQCAH CKFTCHYRCR ALVCLDCCGP RDLGWEPAVE RDTNVDEPVE  
WETPDLSQAE IEQKIKEYNA QINSNLFMSL NKDGSYTGFI KVQLKLVPRV SVPSSKKPPS LQDARRGPGR GTSVRRRTSF YLPKDAVKHL  
HVLSTRRARE VIEALLRKFL VVDDPRKFAL FERAERHGQV YLRKLLDDEQ PLRLRLLAGP SDKALSFVLK ENDSGEVNWD AFSMPPELHNF  
LRILQREEEE HLRQILQKYS YCRQKIQEAL HACPLG

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human RASSF1A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.