

Human CREBZF Recombinant Protein

Catalog Number: PROTQ9NS37

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human CREBZF Recombinant Protein
Description	Human CREBZF Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-354aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol, 1mM DTT, 0.1mM PMSF

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q9NS37, 1-354aa; MGSSHHHHHH SSGLVPRGSH MGSMRHS LTK LLAASGSNSP TRSESPEPAA TC SLPSDLTR AAAGEEETAA
AGSPGRKQQF GDEGELEAGR GSRGGVAVRA PSPEEMEEEA IASLPGEETE DMDFLSGLEL ADLLDPRQPD WHLDPGLSSP GPLSSSGGGS
DSGGLWRGDD DDEAAAAEMQ RFSDLLQRL NGIGGCSSSS DSGSAEKRRR KSPGGGGGGG SGNDNNQAAT KSPRKAAAAA ARLNRLKKKE
YVMGLESVR GLAENQELR AENRELGKRV QALQESRYL RAVLANETGL ARLLSRLSGV GLRLTTS LFR DSPAGDHDYA LPVGKQKQDL
LEEDDSAGGV CLHVDKDKVS VEFCSACARK ASSSLKM

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human CREBZF Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.