

Human OXSM Recombinant Protein

Catalog Number: PROTQ9NWU1

Overview

Product Name	Human OXSM Recombinant Protein
Description	Human OXSM Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 28-459aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9NWU1, 28-459aa; MGSSHHHHHH SGLVPRGSH MGSKRKFFGT VPISRLHRRV VITGIGLVTP LGVGTHLVWD RLIGGESGIV SLVGEEYKSI PCSVAAYVPR GSDEGQFNEQ NFVSKSDIKS MSSPTIMAIG AAELAMKDSG WHPQSEADQV ATGVAIGMGM IPLEVVSETA LNFQTKGYNK VSPFFVPKIL VNMAAGQVSI RYKLGPNHA VSTACTTGAH AVGDSEFRFIA HGDADVMVAG GTDSCISPLS LAGFSRARAL STNSDPKLAC RPFHPKRDGF VMGEGAAVLV LEEYEHAVQR RARIYAEVLG YGLSGDAGHI TAPDPEGEA LRCMAAALKD AGVQPEEISY INAHATSTPL GDAENKAIK HLFKDHAYAL AVSSTKGATG HLLGAAGAVE AAFRTLACY YQKLPPTLNLD CSEPEFDLNY VPLKAQEWKT EKRFGLTNS FGFGGTNATL CIAGL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human OXSM Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.