

Human ADPRHL2 Recombinant Protein

Catalog Number: PROTQ9NX46

Overview

Product Name	Human ADPRHL2 Recombinant Protein
Description	Human ADPRHL2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-363aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 30% glycerol, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9NX46, 1-363aa; MGSSHHHHHH SSGLVPRGSH MGSHEMAAM AAAAGGGAGA ARSLSRFRGC LAGALLGDCV GSFYEAHDTV
DLTSVLRHVQ SLEPDGTPG SERTEALYYT DDTAMARALV QSLLAKEAFD EVDMAHRFAQ EYKKDPDRGY GAGVVTVFKK LLNPKCRDVF
EPARAQFNGK GSYGNGGAMR VAGISLAYSS VQDVQKFARL SAQLTHASSL GYNGAILQAL AVHLALQGES SSEHFLKQLL GHMEDLEGDA
QSVLDARELG MEERPYSRL KKIGELLDQA SVTREEVSE LGNGIAAFES VPTAIYCFLR CMEPDPEIPS AFNSLQRTL IYISLGGDTD TIATMAGAIA
GAYYGMDQVP ESWQQSCEGY EETDILAQSL HRVFQKS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ADPRHL2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.