

Human recombinant FLRT3 protein

Catalog Number: PROTQ9NZU0-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant FLRT3 protein
Description	Human recombinant FLRT3 protein expressed in HEK293 with His-Tag. Sequence domain: 29-528aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q9NZU0, 29-528aa; KSCPSVCRC D AGFIYCND RF LTSIPTGI PE DATTLYLQ NN QINNAGIP SD LKNLLKVER I YLYHNSLDEF PTNLPKYVKE
LHLQENNIRT ITYDSLSKIP YLEELHLD DN SVSAVSIEEG AFRDSNYLRL LFLSRNHLST IPWGLPRTIE ELRLDDNRIS TISSPSLQGL TSLKRLVLDG
NLLNNHGLGD KVFFNLVNL T ELSLVRNSLT AAPVNLPGTN LRKLYLQDNH INRVPPNAFS YLRQLYRLDM SNNNLSNLPQ GIFDDLDNIT
QLILRNNPWY CGCKMKWVRD WLQSLPVKVN VRGLMCQAPE KVRGMAIKDL NAELFDCKDS GIVSTIQITT AIPNTVYPAQ GQWPAPVTQK

PDIKNPCLTK DHQTTGSPSR KTITITVKS SV TSDTIHISWK LALPMTALRL SWLKLGHSPA FGSITETIVT GERSEYLVTA LEPDSPYKVC MVPMETSPLY
LFDETPVCIE TETAPLRMYN PTTTLNREQE KEPYKNPNLP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant FLRT3 protein](#)

For Research Use Only. Not for use in diagnostic procedures.