

Human ABI3 Recombinant Protein

Catalog Number: PROTQ9P2A4

Overview

| | |
|--------------|--|
| Product Name | Human ABI3 Recombinant Protein |
| Description | Human ABI3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-366aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays. |
| Size | Starting from 4ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol |

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q9P2A4, 1-366aa; MGSSHHHHHH SSSLVPRGSH MGSMAELQQL QEFEIPTGRE ALRGNHSALL RVADYCEDNY VQATDKRKAL
EETMAFTTQA LASVAYQVGN LAGHTLRMLD LQGAALRQVE ARVSTLGQMV NMHMEKVARR EIGTLATVQR LPPGQKVIAP ENLPPLTPYC
RRPLNFGCLD DIGHGKDL S TQLSRTG T L S R K S I K A P A T P A S A T L G R P P R I P E P V H L P V V P D G R L S A A S S A S S L A S A G S A E G V G G A P T P K G Q A A P P A P P L
P S S L D P P P P P A A V E V F Q R P P T L E E L S P P P P D E E L P L P L D L P P P P L D G D E L G L P P P P P G F G P D E P S W V P A S Y L E K V V T L Y P Y T S Q K D N E L S F S E G T V I C V
T R R Y S D G W C E G V S S E G T G F F P G N Y V E P S C

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ABI3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.