

Human B3GAT1 Recombinant Protein

Catalog Number: PROTQ9P2W7

Overview

Product Name	Human B3GAT1 Recombinant Protein
Description	Human B3GAT1 Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 28-334aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 0.1M NaCl, 10% glycerol

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9P2W7, 28-334aa; TLAPLLA VHKDEGSDPR RETPPGADPR EYCTSDRDIV EVRTEYVYT RPPPWSDTLP TIHVVTPTYS RPVQKALTR
MANTLLHVPN LHWLVVEDAP RRTPLTARLL RDTGLNYTHL HVETPRNYKL RGDARDPRIP RGTMQRNLAL RWLRETFPRN SSQPGVVYFA
DDDNTYSLEL FEEMRSTRRV SVWPVAFVGG LRYEAPRVNG AGKVVGWKTV FDPHRPFAID MAGFAVNLRL ILQRSQAYFK LRGVKGGYQE
SLLRELVTL NDLEPKAANC TKILVWHTRT EKPVLVNEGK KGFTDPSVEI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human B3GAT1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.