

Mouse TLR2 Recombinant Protein

Catalog Number: PROTQ9QUN7

Overview

Product Name	Mouse TLR2 Recombinant Protein
Description	Mouse TLR2 Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 25-587aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Mouse, Q9QUN7, 25-587aa; ADLQESLSCD ASGVCDGRSR SFTSIPSGLT AAMKSLDLSF NKITYIGHGD LRACANLQVL MLKSSRINTI EGDAFYSLGS LEHLDSLNDH LSSLSSWFG PLSSLKYLNLMGPNYQTLGV TSLFPNLNL QTLRIGNVET FSEIRRIDFA GLTSLNELEI KALSLRNYQS QSLKSIRDIH HLTLLHLSA FLLEIFADIL SSVRYLELRD TNLARFQFSP LPVDEVSSPM KKLAFRGSVL TDESNELLK LLRYILELSE VEFDDCTLNG LGDFNPSESD VVSELGK VET VTIRRLHIPQ FYLFYDLSTV YSLLKVKRI TVENSKVFLV PCSFSQHLKS LEFLDLENL MVEEYLNKNSA CKGAWPSLQT LVLSQNHLSR MQKTGEILLT LKNLTSLDIS RNTFHPMPDS CQWPEKMRFL NLSSTGIRVV KTCIPQTLV LDVSNNNLDS FSLFLPRLQE LYISRNKLT LPDASLFPVL LVMKIRENAV STFSKDQLGS FPKLETLEAG DNHFVCSCEL LSFTMETPAL AQILVDWPDY YLCDSPPRHL GHRLQDARPS VLECHQLEPK SCDKTHTCP CPAPPELLGGP SVFLFPPKPK DTLMISRTPE VTCVVVDVSH EDPEVKFNWY VDGVEVHNAK TKPREEQYNS TYRVVSVLTV LHQDWLNGKE YKCKVSNKAL PAPIEKTISK AKGQPREPQV YTLPPSRDEL TKNQVSLTCL VKGFYPSDIA VEWESNGQPE NNYKTTTPVL DSDGSFFLYS KLTVDKSRWQ QGNVFCSCVM HEALHNHYTQ KSLSLSPGKH HHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse TLR2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.