

Mouse Ak1 Recombinant Protein

Catalog Number: PROTQ9R0Y5

Overview

| | |
|--------------|--|
| Product Name | Mouse Ak1 Recombinant Protein |
| Description | Mouse Ak1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-210aa. Application(s): Enzyme Activity, SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size | Starting from 10ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT |

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, Q9R0Y5, 1-210aa; MGCCVSS EPQEEGGRKT GEKLLKAKII FVVGPGSGK GTQCEKIVQK YGYTHLSTGD LLRAEVSSGS ERGKLSAIM EKGLVPLDT VLDMLRDAML AKVDSSNGFL IDGYPREVKQ GEEFEQKIGQ PTLLEYVDAG AETMTQRLK RGETSGRVDD NEETIKKRLE TYYNATEPVI SFYDKRGIVR KVNAEGTVDT VFSEVCTYLD SLK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Ak1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.