

Human Dkk-3 Recombinant Protein

Catalog Number: PROTQ9UBP4

Overview

Product Name	Human Dkk-3 Recombinant Protein
Description	Human Dkk-3 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 22-350aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9UBP4, 22-350aa; ADPAPAPTAT SAPVKPGPAL SYPQEEATLN EMFREVEELM EDTQHKL RSA VEEMEAEEAA AKASSEVNLA
NLPPSYHNET NTDTKVGNNIT IHVHREIHKI TNNQTGQMFV SETVITSVGD EEGRRSHECI IDEDCGPSMY CQFASFQYTC QPCRGRMLC
TRDSECCGDQ LCVWGHCTKM ATRGSNGTIC DNQRDCQPGL CCAFQRGLLF PVCTPLPVEG ELCHDPASRL LDLITWELEP DGALDRCPA
SGLLCQPHSH SLVYVCKPTF VGSRDQDGEI LLPREVPDEY EVGSFMEEVR QELEDLERSL TEEMALREPA AAAAALLGGE EIHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Dkk-3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.