

Human MORF4L1 Recombinant Protein

Catalog Number: PROTQ9UBU8

Overview

Product Name	Human MORF4L1 Recombinant Protein
Description	Human MORF4L1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-323aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 30% glycerol, 0.15M NaCl, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q9UBU8, 1-323aa; MGSSHHHHHH SGLVPRGSH MGSHPMAPKQD PKPKFQEGER VLCFHGPLYL EAKCVKVAIK DKQVKYFIHY
SGWNKNWDEW VPESRVLKYV DTNLQKQREL QKANQEYAE GKMRGAAPGK KTSGLQQKNV EVKTKKNKQK TPGNGDGGST SETPQPPRKK
RARVDPTVEN EETFMRVEV KVKIPEELKP WLVDWDLIT RQKQLFYLPA KKNVDSILED YANYKKS RGN TDNKEYAVNE VVAGIKEYFN
VMLGTQLLYK FERPYAEIL ADHPDAPMSQ VYGAPHLRL FVRIGAMLAY TPLDEKSLAL LLNYLHDFLK YLAKNSATLF SASDYEVAPP EYHRKAV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human MORF4L1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.