

Human POLL Recombinant Protein

Catalog Number: PROTQ9UGP5

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human POLL Recombinant Protein
Description	Human POLL Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-300aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 30% glycerol, 0.1M NaCl

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9UGP5, 1-300aa; MLMHHQKYLQ RFLGGKREKK QKEACSIPGI GKRMAEKIIE ILESGHLRKL DHISESVPVL ELFSNIWGAG
TKTAQMWWYQQ GFRSLEDIRS QASLTQQAI GLKHYSDFLE RMPREEATEI EQTVQKAAQA FNSGLLCVAC GSYRRGKATC GDVDVLITHP
DGRSHRGIFS RLLDSLREQE FLTDDLVSQE ENGQQQKYLK VCRLPGPGRR HRRLDIIVVP YSEFACALLY FTGSAHFNRS MRALAKTKGM
SLSEHALSTA VVRNTHGCKV GPGRVLPTPT EKDVFRLLGL PYREPAERDW

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human POLL Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.