

## FBXO3 (NM\_012175) Human Recombinant Protein

Catalog Number: PROTQ9UK99

### Overview

Product Name	FBXO3 (NM_012175) Human Recombinant Protein
Description	Recombinant protein of human F-box protein 3 (FBXO3), transcript variant 1
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MAAMETETAPLTLESPTDPLLLILSFLDYRDLINCCYVSRRLSQLSSHDP LWRRHCKKYWLISEEEKTQKNQCWKS L FIDTYS DVGRYIDHYAAIKKAWD  
DLKKYLEPRCP R M V L S L K E G A R E E D L D A V E A Q I G C K L P D D Y R C S Y R I H N G Q K L V P G L L G S M A L S N H Y R S E D L L D V D T A A G G F Q Q R Q G L K Y C L P L T F C I  
HTGLSQYIAVEAAEGRNKNEVFYQCPDQMARNPAAIDMFIIGATFTDWFTSYVKNVVS G G F P I I R D Q I F R Y V H D P E C V A T T G D I T V S V S T S F L P E L S S V H P  
PHYFFTYRIRIEMSKDALPEKACQLDSRYWRITNAKGDVEEVQGPVVG E F P I I S P G R V Y E Y T S C T T F S T T S G Y M E G Y Y T F H F L Y F K D K I F N V A I P R F H M A  
CPTFRVSIARLEMGPDEYEEMEEEEEEEEDEDDSDADMDSEDEDEEERRRRVFDVPIRRRRCSRLF

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



FBXO3 (NM\_012175) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.