

Human Angiotensin-like 2 Recombinant Protein

Catalog Number: PROTQ9UKU9

Overview

Product Name	Human Angiotensin-like 2 Recombinant Protein
Description	Human Angiotensin-like 2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 22-493aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, Q9UKU9, 22-493aa; MGSSHHHHHH SGLVPRGSH MGQEDGFEGT EEGSPREFIY LNRYKRAGES QDKCTYTFIV PQRVTGAIC VNSKEPEVLL ENRVHKQELE LLNNELLKQK RQIETLQQLV EVDGGIVSEV KLLRKESRNM NSRVTQLYMQ LLHEIIRKRD NALELSQLEN RILNQTADML QLASKYKDLE HKYQHLATLA HNQSEIIAQL EEHCQRVPSA RPVPQPPAA PPRVYQPPTY NRIINQISTN EIQSDQNLKV LPPPLPTMPT LTLPSSTDK PSGPWRDCLQ ALEDGHDTS IYLVK PentN RLMQVWCDQR HDPGGWTVIQ RRLDGSVNFF RNWETYKQGF GNIDGEYWLGNENIYWL TNQ GNYKLLV TME DWGRKVF AE YASFRLEPES EYKLR LGRY HGNAGDSFTW HNGKQFTTLD RDHDVYTGNCAHYQKGGWWY NACAHSNLNG VWYRGGHYRS RYQDGVYWAE FRGGSYSLKK VMMIRPNPN TFH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Angiotensin-like 2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.