

Human recombinant ICAM-5 protein

Catalog Number: PROTQ9UMF0

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant ICAM-5 protein
Description	Human recombinant ICAM-5 protein expressed in HEK293 with hIgG-His-Tag. Sequence domain: 23-835aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9UMF0, 23-835aa; LFGLSAVSQE PFWADLQPRV AFVERGGSLW LNCSTNCPRP ERGGLETSLR RNGTQRGLRW LARQLVDIRE
PETQPVCFRR CARRTLQARG LIRTFQRPDR VELMPLPPWQ PVGENFTLSC RVPGAGPRAS LTLTLRGAQ ELIRRSFAGE PPRARGAVLT
ATVLARREDH GANFSCRAEL DLRPHGLGLF ENSSAPREL R TFLSPDAPR LAAPRLLEV G SERPVSC TLD GLFPASEARV YLALGDQNLS
PDVTLEGDAF VATATATASA EQEGARQLVC NVTLGGENRE TRENTIYSF PAPLLTSEP SVSEGMVTV TCAAGAAQALV TLEGVPAAVP

GQPAQLQLNA TENDRRSFF CDATLDVDGE TLIKNRSAEL RVLYAPRLDD SDCPRSWTWP EGPEQTLRCE ARGNPEPSVH CARSDGGAVL
ALGLLGPVTR ALSGTYRCKA ANDQGEAVKD VTLTVEYAPA LDSVGCPERI TWLEGTEASL SCVAHGVPVP DVICVRSGEL GAVIEGLLRV
AREHAGTYRC EATNPRGSAA KNVAVTVEYG PRFEEPSCPS NWTWVEGSGR LFSCEVDGKP QPSVKCVGSG GATEGVLLPL APPDPSRAP
RIPRVLAPGI YVCNATNRHG SVAKTVVSA ESPPEMDEST CPSHQWLEG AEASALACAA RGRPSPGVRC SREGIPWPEQ QRVSREDAGT
YHCVATNAHG TDSRTVTGV EYRPVVAELA ASPPGGVRPG GNFTLTCRAE AWPPAQISWR APPGALNIGL SSNNSTLSVA GAMGSHGGEY
ECAATNAHGR HARRITVRVA GPW

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human recombinant ICAM-5 protein

For Research Use Only. Not for use in diagnostic procedures.