

## Human GITR Ligand/TNFSF18 Recombinant Protein

Catalog Number: PROTQ9UNG2

### Overview

|              |   |
|--------------|---|
| Product Name | Human GITR Ligand/TNFSF18 Recombinant Protein   |
| Description  | Human GITR Ligand/TNFSF18 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 72-199aa. Application(s): SDS-PAGE. |
| Size         | Starting from 20ug  |
| Tag          | Non-Tagged  |
| Form         | Liquid  |
| Source       | E. coli   |
| Formulation  | 10mM Sodium Citrate buffer (pH 3.5) containing 1mM DTT, 10% glycerol  |

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q9UNG2, 72-199aa; MQLETAKEPC MAKFGPLPSK WQMASSEPPC VNKVSDWKLE ILQNGLYLIY GQVAPNANYN DVAPFEVRLY  
KNKDMIQTLT NKSQIQNVGG TYELHVGDTI DLIFNSEHQV LKNNTYWGII LLANPQFIS

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human GITR Ligand/TNFSF18 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.