

## SAPS1 (PPP6R1) (NM\_014931) Human Recombinant Protein

Catalog Number: PROTQ9UPN7

### Overview

|              |   |
|--------------|---|
| Product Name | SAPS1 (PPP6R1) (NM_014931) Human Recombinant Protein              |
| Description  | Recombinant protein of human SAPS domain family, member 1 (SAPS1) |
| Size         | 20 µg   |
| Tag          | C-Myc/DDK   |
| Form         | Frozen Solution in PBS Buffer                                     |
| Source       | HEK293T   |
| Formulation  | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol              |

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MFWKFDLHTSSHLDTLEREDLSPELLDEEDVLQECKVVRKLLDFLLQPPHLQAMVAWVTQEPPDSSGEERLRYKYPVACEILTSQVDPQINDALGADE  
SLLNRLYGFLQSTGSLNPLLASFFSKVMGILINRKTDLVSLFRKKDDFVDFLLQHGTSAIMDLLLRLLTCVERPQLRQDVVNWLNEEKIVQRLIEQIHPSK  
DENQHSNASQSLCDIIRLSREQMIQVQDSPEPDQLLATLEKQETIEQLLSNMFEQEQQSVIVSGIQVLLTLEPRRPRSESVTVNSFFSSVDGQLELLAQQ  
ALESTVSSVGHALHALRPLRSCFHQLLLEPPKLEPLQMTWGMPLAPPLGNTRLHVVKLLASALSANDAALTHELLALDVPNTMLDLFFHYVFNFLHAQVEG  
CVSTMLSLGPPPDSSPETPIQNPVVKHLLQQCRLVERILTSWEENDRVQCAGGPRKGYMGHLTRVAGALVQNTTEKGNPNAEQLRQLLKELPSEQEQWE  
AFVSGPLAETNKKNMVDLVNTHHLHSSDDEDDRLKEFNFPPEAVLQAFMDFQMQRMTSAFIDHFGFNDEEFGEQEESVNAPFDKTANITFSLNAD  
DENPNANLLEICYKDRIQQFDDDEEEDEEEAQQSGESDGEDGAWQGSQQLARGARLQPPGVRSGGSTDSEDEEEDEEEDEEEDGIGCAARGGATPL  
SYSPGPQPPGPSWTATFDVPVPTDAPTSRVS GEEELHTGPPAPQGPLSVPQGLPTQSLASPPARDALQLRSQDPTPPSAPQEATEGSKVTEPSAPCQAL  
VSGIDLQATFHGIRSAPSSSDSATRDPSTVSPASGAHQPPQTTEGEKSPEPLGLPQSQAQALTPPPIPNGSAPEGPASPGSQ

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



SAPS1 (PPP6R1) (NM\_014931) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.