

## Human PA2G4 Recombinant Protein

Catalog Number: PROTQ9UQ80

### Overview

Product Name	Human PA2G4 Recombinant Protein
Description	Human PA2G4 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-394aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, Q9UQ80, 1-394aa; MSGEDEQQEQ TIAEDLVVTK YKMGGDIANR VLRLSVEASS SGVSVLSLCE KGDAMIMEET GKIFKKEKEM  
KKGIAFP TSI SVNNCVCHFS PLKSDQDYIL KEGDLVKIDL GVHVDGFIAN VAHTFVVDVA QGTQVTGRKA DVIKAAHLCA EAALRLVKPG  
NQNTQVTEAW NKVAHSFNCT PIEGMLSHQL KQHVIDGEKT IIQNPTDQK KDHEKAEFEV HEVYAVDVLV SSGEGKAKDA GQRTTIYKRD  
PSKQYGLKMK TSRAFFSEVE RRFDAMPFTL RAFEDEKKAR MGVVECAKHE LLQPFNVLYE KEGEFVAQFK FTVLLMPNGP MRITSGPFEP  
DLYKSEMEVQ DAELKALLQS SASRKTQKKK KKKASKTAEN ATSGETLEEN EAGDLEHHHH HH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human PA2G4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.