

## Mouse Ak2 Recombinant Protein

Catalog Number: PROTQ9WTP6

### Overview

Product Name	Mouse Ak2 Recombinant Protein
Description	Mouse Ak2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-239aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 5ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.5) containing 10% glycerol, 1mM DTT

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q9WTP6, 1-239aa; MGSSHHHHHH SSGLVPRGSH MGSHMAPNVL ASEPEIPKGI RAVLLGPPGA GKGQTQAPKLA ENFCVCHLAT  
GDMLRAMVAS GSELGKLLKA TMDAGKLVSD EMVVELIEKN LETPCKNGF LLDGFPRTVR QAEMLDDLME KRKEKLD SVI EFSIQD SLLI  
RRITGRLIHP KSGRSYHEEF NPPKEPMKDD ITGEPLIRRS DDNEKALKTR LEAYHTQTTP LVEYYRKRGI HCAIDASQTP DIVFASILAA FSKATCKDLV  
MFI

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Ak2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.