

Human CCL26/Eotaxin-3 Recombinant Protein

Catalog Number: PROTQ9Y258

Overview

Product Name	Human CCL26/Eotaxin-3 Recombinant Protein
Description	Human CCL26/Eotaxin-3 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 24-94aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9Y258, 24-94aa; MTRGSDISK TCCFQYSHKPL PWTWVRSYEF TSNCSQRAV IFTTKRGKKV CTHPRKKWVQ KYISLLKTPK QL

[?/description_after_attributes?](#)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CCL26/Eotaxin-3 Recombinant Protein