

B3GALT5 (NM_033173) Human Recombinant Protein

Catalog Number: PROTQ9Y2C3

Overview

Product Name	B3GALT5 (NM_033173) Human Recombinant Protein
Description	Purified recombinant protein of Homo sapiens UDP-Gal:betaGlcNAc beta 1,3-galactosyltransferase, polypeptide 5 (B3GALT5), transcript variant 5
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MAFPKMRLMYICLLVLGALCLYFSMYSLNPFKEQSFVYKKGDNFLKLPDTCRQTPPFLVLLVTSSHKQLAERMAIRQTWGKERTVKGKQLKTFLLGT
TSSAAETKEVDQESQRHGDIIQKDFLDVYYNLTLLKMMGIEWVHRFCPQA AFVMKTDSDMFINVDYLT ELLKKNRTRRFFTGFLKLNFPPIRQPFSKW
FVSKSEYPWDRYPPFCSGTGYVFGDVASQVYNVSKSVPIYKLEDVVFVGLCLERLNIRLEELHSQPTFFPGGLRFSVCLFRRIVACHFIKPRLLDYWQALE
NSRGEDCPPV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



B3GALT5 (NM_033173) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.