

Human recombinant CARHSP1 protein

Catalog Number: PROTQ9Y2V2-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant CARHSP1 protein
Description	Human recombinant CARHSP1 protein expressed in E.coli with His-Tag. Sequence domain: 1-147aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E.coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.1M NaCl

Concentration

0.5mg/ml (determined by BCA assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9Y2V2, 1-147aa; MGSSHHHHHH SSGLVPRGSH MSSEPPPPPPQ PPTHQASVGL LDTPRSRERS PSPLRGNVVP SPLPTRRTRT
FSATVRASQG PVYKGVCKCF CRSKGHGFI PADGGPDIFL HISDVEGEYV PVEGDEVYK MCSIPPKNEK LQAVEVWITH LAPGTKHETW SGHVISS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant CARHSP1 protein](#)

For Research Use Only. Not for use in diagnostic procedures.